Simple Simon UA Plus [ZAF-702] Urine test strip analyzer

- Fast and quiet
- 60 test / hour
- Barcode scanner
- Thermal receipt printer
- Ethernet and Wifi
- 7" touch display
- Speaker
- USB and Bluetooth
- Poct1A Generic access point
- Continua FHIR Gateway
- CE

USD 1750

(VAT will be added)

The urine strip analyzer, Simple Simon UA Plus from Zafena is a flexible, modular, analyzer that consists of three main parts; The urine strip reader (Meter), The POC-Workstation (IT-support) and a barcode scanner (Scanner).

The POC-Workstation hosts up to three simultaneous cable connected and multiple BT connected meters. It is very easy to connect new meters through the patented instrument-tag that connect and pair the devices without extended IT knowledge or issues.

The meter has a sled for placing the urine strip, one ON/OFF button and one START button. The meter is connected to the POC-Workstation through a USB cable and can be placed at a distance from the POC-Workstation or hung on the wall. The interaction with the analyzer is through the 7" large color touch display or the attached barcode scanner.

The built in receipt printer is for printing results, but foremost to print out lab-orders, including laboratory-ID and requested analytes.

- *Visual verification*
The 7" display present the results, but also the actual detected color changes for the operator to visually confirm. The colors of a dry strip is compared to the colors detected by the analyzer.
- *Step by step instructions*
The step by step instructions shown on the display will lead the operator through the measurement
- *No settings menu*
Settings are changed through scanning of barcodes. No worries to accidentally change vital settings.
- *Arbitrary or SI units*
The presentation of the results can be changed between SI and arbitrary units.
- *Multi-connect*
Connect more than one meter to the POC-Workstation to increase production. Or widen the application by adding more suitable analytes to the same workstation.
- *Printer*
Print out of results and lab orders direct from EHR. IT also supports ipp printing to local network printers,

CYBOW 8A

[ZAF-130-8] 8 parameter urine test strip

- 100 test strips per container
- Humidity protection in the container
- 2 year shelf life from manufacturing
- 6 month once opened
- Ascorbic acid detection
- CE

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Product number	Description
ZAF-702	Simple Simon UA Plus Urine strip analyzer. Including ZAF-552 POC-Workstation and ZAF-130-1 Meter
ZAF-130-1	Simple Simon UA - Urine strip meter
ZAF-130-8	CYBOW 8A - Eight parameter urine strip
ZAF-552	POC-Workstation for connectivity and interaction. Including ZAF-502 Workstation, ZAF-504 Barcode reader, ZAF-130-12 Power adaptor, ZAF-509 Network cable, ZAF-130-13 thermal printer paper roll, ZAF-505-2 Micro SD card
ZAF-130-11	Sled for Simple Simon UA
ZAF-130-13	Thermal printer paper roll (5 pcs)