

Boditech Control

Boditech Control	Catalog No.	Boditech Control	Catalog No.
ichroma CRP control	CFPO-2	Boditech Anti-HBs Control	CFPO-144
ichroma MAU control	CFPO-4	Boditech AFP Control	CFPO-248
ichroma iFOB control L	CFPO-14-1	Boditech CEA Control	CFPO-246
ichroma iFOB control H	CFPO-14-2	Boditech CK-MB Control	CFPO-243
hemochroma Plus control	CFPO-29	Boditech Cortisol Control	CFPO-236
ichroma vitamin-D control	CFPO-79	Boditech FSH Control	CFPO-230
Boditech Tumor marker Control	CFPO-94	Boditech HCG Control	CFPO-232
Boditech Hormone Control	CFPO-95	Boditech LH Control	CFPO-234
Boditech HbA1c Control	CFPO-96	Boditech Myoglobin Control	CFPO-244
Boditech PCT Control	CFPO-97	Boditech NTproBNP Control	CFPO-245
Boditech Cardiac Control	CFPO-98	Boditech PRL Control	CFPO-226
Boditech Ferritin Control	CFPO-99	Boditech Progesterone Control	CFPO-238
Boditech CRP Control	CFPO-100	Boditech T3 Control	CFPO-240
Boditech D-Dimer Control	CFPO-101	Boditech T4 Control	CFPO-237
Boditech RF IgM Control	CFPO-103	Boditech Testosterone Control	CFPO-239
Boditech Anti-CCP Control	CFPO-125	Boditech Tn-I Control	CFPO-241
Boditech HBsAg Control	CFPO-142	Boditech Tn-I Plus Control	CFPO-212
Boditech Anti-HCV Control	CFPO-143	Boditech TSH Control	CFPO-228

Boditech Calibrator

Boditech Calibrator	Catalog No.	Boditech Calibrator	Catalog No.
Boditech Tumor marker Calibrator	CFPO-2	Boditech FSH Calibrator	CFPO-256
Boditech Hormone Calibrator	CFPO-4	Boditech HCG Calibrator	CFPO-258
Boditech HbA1c Calibrator	CFPO-14-1	Boditech LH Calibrator	CFPO-260
Boditech PCT Calibrator	CFPO-14-2	Boditech Cortisol Calibrator	CFPO-262
Boditech Cardiac Calibrator	CFPO-29	Boditech T4 Calibrator	CFPO-263
Boditech Ferritin Calibrator	CFPO-79	Boditech Testosterone Calibrator	CFPO-265
Boditech CRP Calibrator	CFPO-94	Boditech T3 Calibrator	CFPO-266
Boditech D-Dimer Calibrator	CFPO-95	Biotech CK-MB Calibrator	CFPO-269
Biotech AFP Calibrator	CFPO-274	Biotech Myoglobin Calibrator	CFPO-270
Boditech PRL Calibrator	CFPO-252	Biotech NTproBNP Calibrator	CFPO-271
Boditech TSH Calibrator	CFPO-254	Biotech CEA Calibrator	CFPO-272

AFIAS-6

Versatile desktop immuno-analyzer

- Automated desktop analyzer
- 10 μ L or 30 μ L finger prick blood
- 6 different tests on 6 different channels
- Bay random access (2 bays)

All-in-one cartridge format
Just load sample & press start

AFIAS

AFIAS-1

Doctor's reliable associate

- Automated desktop analyzer
- LIS / HIS connectivity
- All-in-one cartridge
- 10 μ L or 30 μ L finger prick blood
- Walk away function

Specifications

QC / Calibration

2 calibrators & 2 controls
Internal & external controls

Interface

7" touch screen
Built-in thermal printer
RS 232 / USB port / Ethernet / SD card slot
LIS / HIS compatible

Memory

5,000 patient results
5,000 QC results
500 ID chips
100 User IDs

Throughput

Up to 36 tests / hour

Analyzer

420 x 336 x 293 mm
15.1 kgs
100 – 240V AC, 50-60Hz
Internal temperature control

Specifications

QC / Calibration

2 calibrators & 2 controls
Internal & external controls available

Interface

5" touch screen
Built-in thermal printer
RS 232 / USB port / Ethernet / SD card slot
LIS / HIS compatible

Memory

5,000 patient results
5,000 QC results
500 ID chips
100 User IDs

Throughput

Up to 10 tests / hour

Analyzer

320 x 204 x 180 mm
3.9 kgs
100 – 240V AC, 50-60Hz
Internal temperature control