

Godkänt av:	Version: 2019-06-11	
Datum:	Författare:FW	Sidan 1 av 2

## Available instrument-tags

### Description

The document describes all instrument-tag bar codes.


### Instrument-tags

An instrument-tag consists of two or three type of data:  
 "LMC" indicates that an instrument-view should be opened.  
 "-XXX1" tell the Plus-Screen which view to start.  
 "-0001" (optional) customer specific meter identification.


#### Simple Simon PT




#### Simple Simon UA




#### TD-3128B (Taidoc)




#### TD-1107B (Taidoc)




#### TD-4289 (Taidoc)




#### TD-2555B (Taidoc)




#### TD-1241B (Taidoc)



#### TD-8255 (Taidoc)



#### BYH-03 (Beryl scales)



#### Diaspect Hemoglobin



#### Diaspect TM



#### HemoControl



#### Elitech MicroSed R




#### Vacurette SRT/10



#### Terumo V10



#### LumiraTek C



#### QuikReadGo





#### DCA Vantage




Godkänt av:	Version: 2019-06-11	
Datum:	Författare:FW	Sidan 2 av 2

**Clinitek Status**  
(Urine strip analyzer)  
  
LMC-BCS1


**HemoQue WBC Diff**  
(WBC)  
  
LMC-HWD1


**HemoQue Glucose 201+**  
(B-glucose)  
  
LMC-HGL1


**HemoQue Hb 201+**  
(Hb)  
  
LMC-HHB1


**Boule Exigo H400**  
(Hematology)  
  
LMC-BEX1

**MicroINR**  
(PT/INR Quick)  
  
LMC-IMI1

**Boule C200**  
(Hematology)  
  
LMC-BEX2


**Uryxxon Relax**  
(Urine strip analyzer)  
  
LMC-MUR1


**I-Chroma**  
(CRP)  
  
LMC-BIC1


**Protime InRhythm**  
(PT/INR Quick)  
  
LMC-IIR1

**Simple Simon AP**  
(APTT)  
  
LMC-ZAP1


**Simple Simon DD**  
(D-dimer)  
  
LMC-ZDD1


**HemCheck Helge 1**  
(Hemolysis)  
  
LMC-HCH1

**Affinion 2**  
(HbA1C, Lipid, CRP)  
  
LMC-AAF2

**Quo-test (EKF)**  
(HbA1C)  
  
LMC-EQT1

**LDX Cholestech**  
(Lipid panel)

**Peak Flow Meter**  
(PEF, FEV1)  
  
**Bluetooth**  
(The unit MAC-adress)

**EVOLV (Omron)**  
(Blood pressure)  
  
**Bluetooth**  
(The unit MAC-adress)

**The Cube (Eurolyser)**  
(In progress)

**M2 (Omron)**  
(Blood pressure)

**U18 (Omron)**  
(Blood pressure)

**RS7 Intelli IT (Omron)**  
(Blood pressure)

**M10IT (Omron)**  
(Blood pressure)

**AccuCheck (Roche)**  
(glucose)

**MIT Elite Plus (Omron)**  
(Blood pressure)

**Xpresia Stride**  
(PT/INR)

**Micropoint (qLabs)**  
(PT/INR)

**TD-4216B (Taidoc)**  
(glucose, ketone, lactate,  
Cholesterol, Uric acid)

**TD-1241 (Taidoc)**  
(non touch temp)