

MicroINR Instructions

Approved by:	Version: 20230828	
Date:	Author: FW	Page 1 by 1

iLine MicroINR Instructions

The document describes the connection and labeling of the meter, to work with the Zafena POC-Workstation. Physical connection is between USB-A and USB-Mini-B.

Produkt kod kabel: **ZAF-520-4**



Instrument-tag: LMC-IMI1



Print and attach the instrument-tag on to the meter. In some cases, to enhance the work flow, it might work better to attach the tag in close proximity to the meter.

Communication settings

There are no communication settings to set on the meter. When the measurement is complete and the result is shown on the meter's display, it is also sent digitally via the cable to the POC-Workstation.