

Quo-Test®



Quo-Test®

Easy and reliable HbA1c measurement

- ◆ Fully automated analyser
- ◆ Fast and accurate
- ◆ Easy to use
- ◆ Results in just 4 minutes

Simple, accurate and reliable HbA1c measurement

<p>Easy to use</p>	<ul style="list-style-type: none"> • Intuitive, step by step on screen instructions minimise training time • Step by step instructions on display • Multi-lingual menu • Just 4µl blood required from finger prick or venous sample
<p>Efficient data handling</p>	<ul style="list-style-type: none"> • Barcode reader included. Only one scan per lot of cartridges required • Stores up to 7000 results for efficient patient data management • User selectable dual reporting. User can select % DCCT, IFCC mmol/mol, eAG mg/dl or eAG mmol/l. • USB Port • Optional label printer available
<p>Fast and accurate</p>	<ul style="list-style-type: none"> • Results within 4 minutes • Boronate affinity method which is recognized as interference free • Measuring range: 4–15% A1c DCCT • Imprecision: CV <3% at 7% A1c DCCT • Unaffected by Hb variants, which do not result in reduced erythrocyte life span • Traceable to the IFCC reference method
<p>Ergonomic design</p>	<ul style="list-style-type: none"> • Compact: 205 mm(h) x 205 mm(w) x 135 mm(d) • Lightweight: 1.3kg

Catalog No.	Name
7045-0011-0070	• Quo-Test Analyzer Pack



1 Collect 4 µl blood sample.



2 Insert sample collector into the Quo-Test cartridge.



3 Insert cartridge into analyser. Results displayed within 4 minutes.

Distributed by



Quotient Diagnostics:
ISO13485: 2012/2003 Certified.
Registration ref: MD 505902

Manufacturer
Quotient Diagnostics Ltd.
Russell House, Molesey Road
Walton-on-Thames, Surrey, UK,
KT12 3PE

www.ekfdiagnostics.com

EKF | Diagnostics for life